



Soft Tissue-Supported Full-Arch Guide Rx

Please fill out the form completely
Contact us at **1-800-216-5407** with any questions
Please allow two weeks from case approval for delivery

General Info

Doctor Name:

Business Name and Address:

Patient First Initial and Last Name:

Best Email:

Best Phone Number:

What is included in soft tissue supported guide cases?

- Treatment planning assistance provided by a DENTIST who has clinical knowledge of these procedures.
- Smile design simulation
- Drill report and video walk-through of the case
- 3D printed soft tissue supported pin-retained drill guide (includes backup duplicate)
- 3D printed models of patient to verify fit and to perform mock surgery if desired

Printed Guide per Arch:

\$850

Digital Design Only:

\$500

*****DISCLAIMER*****

These are complex procedures and the treating doctor must fully understand how to perform and manage any and all aspects of the treatment, including complications. Transcend Digital Dental Lab does not warrant, guarantee, or assume responsibility for the accuracy of your surgery or the final clinical outcome. We exist as a resource for treatment planning assistance, guide design, and fabrication and, we are happy to review any and all aspects of the case for your benefit. However, YOU as the treating doctor are the one responsible for the clinical results.

Please sign here that you have read, understood, and agreed to the disclaimer above:

Printed name:_____ Signature_____ Date:_____

Case Information

Implant System:

Surgical Kit:

Circle jaw: Upper jaw Lower Jaw Both

Number of Implants per Arch: 4 5 6 7 8

Can we plan with Vertical Sinus Lifts(bumps, tents)? Yes No

Implant Positioning: Straight Tilted No Preference

Fixation Pins for Guide: Yes No

(We recommend using Blue Sky Bio Pins)

If using another type of pin, please indicate brand and dimensions of the pin

Metal sleeves for guide? (\$20 per sleeve) Yes No

Metal sleeves not recommended when using a fully-guided kit or a keyed guided system

Metal Sleeves recommended if any cutting portion of the drill will come in contact with the guide

Immediate Load Case? Yes No

Include printed immediate load prosthesis? Yes No

(\$300 per arch w/ duplicate; \$500 per arch if milled)

Restoration designed with premade 5mm holes for the emergence of temp cylinders.

Intaglio is designed either flat or ovate for expedited conversion. Printed permanent crown and bridge resin used unless otherwise indicated.

Backup 3D printed denture (Add \$200 per arch) Yes No

(In case adequate primary stability is not achieved)

Shade: B1 A2 A3

Permission to Document and Share case as Lecture Content? Yes No

(We will maintain HIPAA compliance and hide pt. identity)

Checklist of Data to Send

Scan 1: CBCT scan of patient wearing ideal denture with radiographic markers.

- Please position the patient **IN OCCLUSION AT THE DESIRED VDO** using a **CHIN REST** rather than bite fork.
- **Scan 2:** CBCT of the denture by itself with the radiographic markers (dual scan)
- Pictures:
 - Full face at repose
 - Full face with highest smile
 - Full face with retractors in place
 - Full face from a lateral view
- ***Please ensure all digital models are in proper occlusion
OR
- Physical models/impressions come with bite registration

Additional Notes Regarding Case:

Please upload digital data to **BioBigBox.com** via email **transcenddigitallab@gmail.com**

Any physical models can be sent to:

Transcend Digital Dental Lab
4949 Century Street NW
Huntsville, AL 35816